



# The Scarban® assortment


## LIGHT



## ELASTIC




## ELASTIC (special products)



Size definition: contour palm of the hand

S	10 - 11 cm	S	18 - 19 cm
M	12 - 13 cm	M	20 - 21 cm
L	14 - 15 cm	L	22 - 23 cm
XL	16 - 17 cm	XL	24 - ++ cm



## Indicated for:

- Prevention and treatment of hypertrophic and keloid scars in children
- Treatment of burns (as monotherapy or in combination with pressure).

## Scarban® effectiveness:

**Case I: Burn scar**  
Patient: kid  
Cause: burn

Before



During



After



**Case II: Burn scar**  
Patient: male  
Cause: car crash  
(In combination with pressure)

Before



After




## Product + brand features:

- Easy application + usage
- For the silicone sheets: washable, reusable, cost efficient, double elastic + UV protective (UPF50)
- Excellent semi-occlusiveness
- Extended range available (different types of silicone, different shapes + forms)
- Long lasting treatment per package (depending size & maintenance)



## ELASTIC SILICONE SHEETS



Double  
Elastic



Thickness  
1,4  
mm



Re-Usable

### Where:

- Large surfaces
- Hypertrophic/keloid scars
- Large burn scars
- Hands

## LIGHT SILICONE SHEETS



Double  
Elastic



Thickness  
0,4  
mm



Re-Usable

### Where:

- Smaller surfaces
- Postoperative/linear incisions
- Face
- Fingers